

BURGUNDY CHROME PLATED DUAL TUBE BINAURALS

Réf. 507923



BURGUNDY CHROME PLATED DUAL TUBE BINAURALS

Réf. 507923

FEATURES			
Medical device	Yes	Classe I	
Accuracy measurement			Classe I
Marquage		Barcode	3700446009601
Guarantee (month):	0	Customs Code	90181990
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,080
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:47 - This article is available on our website : www.gsh-med.fr with the following reference 507923

All characteristics are given for information only and may be changed