

FLOATING DIAPHRAGM/BLACK RETAINING RING CARDIO SONIC

Réf. 507934



FLOATING DIAPHRAGM/BLACK RETAINING RING CARDIO SONIC

Réf. 507934

Spare part for Cardio Sonic® stethoscopes

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446009656
Guarantee (month):	0		Customs Code	90181990
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :			Width packaging (cm):	
Length (cm) :			Length packaging (cm) :	
Height (cm) :			Height packaging (cm) :	
Weight (kg) :			Weight packaging (kg) :	0,005
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



31/03/2022 01:47 - This article is available on our website : www.gsh-med.fr with the following reference 507934

All characteristics are given for information only and may be changed