

DIAPHRAGM/BLACK RETAINING RING FOR LAUBRY

Réf. 507936



DIAPHRAGM/BLACK RETAINING RING FOR LAUBRY

Réf. 507936

Spare part for Laubry® stethoscopes

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446009670
Guarantee (month):	0	Customs Code	90181990
Sterile	No		
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	4,9
Length (cm) :		Length packaging (cm) :	4,9
Height (cm) :		Height packaging (cm) :	0,4
Weight (kg) :		Weight packaging (kg) :	0,005
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:47 - This article is available on our website : www.gsh-med.fr with the following reference 507936

All characteristics are given for information only and may be changed