

DIAPHRAGME/GREY RETANNING RING FOR PULSE & DUAL PUSLE



Réf. 507937



DIAPHRAGME/GREY RETANNING RING FOR PULSE & DUAL PUSLE



Réf. 507937

Spare part for Pulse® & Dual Pulse® stethoscopes

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446009687
Guarantee (month):	0		Customs Code	90181990
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :			Width packaging (cm):	
Length (cm) :			Length packaging (cm) :	
Height (cm) :			Height packaging (cm) :	
Weight (kg) :			Weight packaging (kg) :	0,005
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



15/11/2022 12:46 - This article is available on our website : www.gsh-med.fr with the following reference 507937

All characteristics are given for information only and may be changed