

PEDIATRE BELL BLACK RETAINING RING CARDIO PRESTIGE/MAGISTER FOR CHILD

Réf. 507943



PEDIATRE BELL BLACK RETAINING RING CARDIO PRESTIGE/MAGISTER FOR CHILD

Réf. **507943**

Spare part for Cardio Prestige® stethoscopes

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446009724
Guarantee (month):	0		Customs Code	90181990
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :	0,005	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



31/03/2022 01:48 - This article is available on our website : www.gsh-med.fr with the following reference 507943

All characteristics are given for information only and may be changed