

NON CHILL RING (BLUE) FOR BELL OF DUAL PULSE II SPENGLER



Réf. 507944

NON CHILL RING (BLUE) FOR BELL OF DUAL PULSE II SPENGLER



Réf. 507944

FEATURES			
Medical device	No		
Marquage		Barcode	3700446030087
Guarantee (month):	0	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,001
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:48 - This article is available on our website : www.gsh-med.fr with the following reference 507944

All characteristics are given for information only and may be changed