

## SINGLE TUBBING BLACK WITH BINAURAL ASSEMBLY FOR PEDIASONC SPENGLER

Réf. 507950



# SINGLE TUBBING BLACK WITH BINAURAL ASSEMBLY FOR PEDIASONC SPENGLER

Réf. 507950

Spare part for Pédiasonic® stethoscopes

FEATURES			
Medical device	Yes	Classe I	
Marquage		Accuracy measurement	Classe I
Guarantee (month):	0	Barcode	3700446011277
Sterile	No	Customs Code	90181990
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,070
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:48 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 507950

All characteristics are given for information only and may be changed