

LIAN NM CLINIC / LARGE ADULT (L)

Réf. 520660



LIAN NM CLINIC / LARGE ADULT (L)



Réf. 520660

High resistance model with Clinic® cuff, quick coupling and clamp for liquid decontamination. Robust, chromium-plated metal casing with foam manometer cover. This model is especially appreciated by emergency services. Supplied in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446024826
Guarantee (month):	60	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,662
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	85	Humidity max	85



31/03/2022 02:30 - This article is available on our website : www.gsh-med.fr with the following reference 520660

All characteristics are given for information only and may be changed