

AUTOTENSIO BPM ARM CUFF ADULT (M/L) BLUE SPENGLER

Réf. SPG440-



AUTOTENSIO BPM ARM CUFF ADULT (M/L) BLUE SPENGLER



Réf. **SPG440-**

AUTOTENSIO® upper arm model. For occasional self-monitoring. - LCD screen, - W.H.O risk indicator - Irregular heart beat detection, - Low batteries indicator, - 2 x 60 memories with date and time, - Purple cuff - Upper arm circumference 22-36 cm size M. Supplied with 4 LR06 AA batteries and a storage bag. Optional power adapter.

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0459	Barcode	3700446022501
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	13,0
Length (cm) :		Length packaging (cm) :	17,0
Height (cm) :		Height packaging (cm) :	10,0
Weight (kg) :		Weight packaging (kg) :	0,594
SPECIFIC STORAGE CONDITION			
Temperature min		Temperature max	55
Humidity min		Humidity max	90



08/12/2022 12:39 - This article is available on our website : www.gsh-med.fr with the following reference SPG440-

All characteristics are given for information only and may be changed