

SUNTECH CT50 BP STANDARD SPO2 PRINTER

Réf. DSU222202

SunTech Medical®

GARANTIE
2
ANS
YEARS
GUARANTEE

CE
0413

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE II



SUNTECH CT50 BP STANDARD SPO2 PRINTER

Réf. DSU222202

The SUNTECH CT 50 monitor is a condensed version of the latest technology.
 The world reference for measuring vital signs.
 Nellcor SpO2 technology, optional HIS, USPx4, Ethernet LAN and WIFI connectivity.
 3 blood pressure measurement modes: Automatic, Auscultatory and average mode.
 8-inch colour touch screen, optional barcode reader, optional printer and optional wheeled base

FEATURES					
Medical device	Yes	Classe IIa	Accuracy measurement	Classe II	
Marquage	CE0413		Barcode	3700446025120	
Guarantee (month):	24		Customs Code		
Sterile	No				
INFORMATIONS		Main unit	Packaged unit		
Width (cm) :		Width packaging (cm) :			
Length (cm) :		Length packaging (cm) :			
Height (cm) :		Height packaging (cm) :			
Weight (kg) :		Weight packaging (kg) :		5,020	
SPECIFIC STORAGE CONDITION					
Temperature min		(N/A)	Temperature max		(N/A)
Humidity min		(N/A)	Humidity max		(N/A)



31/03/2022 01:56 - This article is available on our website : www.gsh-med.fr with the following reference DSU222202

All characteristics are given for information only and may be changed