

DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM)

Réf. SPG440P



DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM)



Réf. **SPG440P**

All characteristics are given for information only and may be changed

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0459	Barcode	3661040214207
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	3,730
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:37 - This article is available on our website : www.gsh-med.fr with the following reference SPG440P

DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM)



Réf. SPG440P

COMPOSITION

Code	Descriptive	Quantity
SPG440-	Autotensio BPM arm cuff Adult (M/L) Blue Spengler	6