

EXAMINATION LAMP LED 3W HOLTEX - BLACK



Réf. OSLAMP01



EXAMINATION LAMP LED 3W HOLTEX - BLACK



Réf. OSLAMP01

FEATURES			
Medical device	No		
Marquage		Barcode	3661040212814
Guarantee (month):	12	Customs Code	94054090000
Sterile	No		
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	12,0
Length (cm) :		Length packaging (cm) :	33,0
Height (cm) :		Height packaging (cm) :	111,0
Weight (kg) :		Weight packaging (kg) :	7,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/08/2022 21:09 - This article is available on our website : www.gsh-med.fr with the following reference OSLAMP01

All characteristics are given for information only and may be changed