

EXAMINATION LAMP LED Q3 WHITE (WITH WHITE METAL STAND 5 WHEELS)

Réf. OSLAMP02



EXAMINATION LAMP LED Q3 WHITE (WITH WHITE METAL STAND 5 WHEELS)



Réf. OSLAMP02

FEATURES			
Medical device	No		
Marquage		Barcode	3661040212821
Guarantee (month):	12	Customs Code	94054090000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	12,0
Length (cm) :		Length packaging (cm) :	33,0
Height (cm) :		Height packaging (cm) :	111,0
Weight (kg) :		Weight packaging (kg) :	7,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:47 - This article is available on our website : www.gsh-med.fr with the following reference OSLAMP02