

## REPLACEMENT SPOON FOR MOBI BLUEBERRY

Réf. 513721



# REPLACEMENT SPOON FOR MOBI BLUEBERRY



Réf. 513721

FEATURES			
Medical device	No		
Marquage		Barcode	3700446031473
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	2,5
Length (cm) :		Length packaging (cm) :	7,0
Height (cm) :		Height packaging (cm) :	2,5
Weight (kg) :		Weight packaging (kg) :	0,006
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:53 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 513721

All characteristics are given for information only and may be changed