

PRINTER PAPER CT50

Réf. DSU222209

PRINTER PAPER CT50

Réf. DSU222209

FEATURES			
Medical device	No		
Marquage		Barcode	3661040214276
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	5,0
Length (cm) :		Length packaging (cm) :	20,0
Height (cm) :		Height packaging (cm) :	5,0
Weight (kg) :		Weight packaging (kg) :	0,188
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



06/04/2022 12:34 - This article is available on our website : www.gsh-med.fr with the following reference DSU222209