

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 CORAL

Réf. SPG345P



AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 CORAL



Réf. **SPG345P**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0459	Barcode	3700446041601
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:00 - This article is available on our website : www.gsh-med.fr with the following reference SPG345P

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 CORAL



Réf. SPG345P

COMPOSITION

Code	Descriptive	Quantity
SPG345	AUTOTENSIO Wrist - Wrist Blood Pressure Monitor Spengler Coral	6