

AUTOTENSIO UPPER ARM (M/L) BLOOD PRESSURE MONITOR SPENGLER CORAL

Réf. SPG445



CE
0459



AUTOTENSIO UPPER ARM (M/L) BLOOD PRESSURE MONITOR SPENGLER CORAL



Réf. **SPG445**

FEATURES			
Medical device	No		
Marquage	CE0459	Barcode	3700446041571
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	13,0
Length (cm) :		Length packaging (cm) :	17,0
Height (cm) :		Height packaging (cm) :	10,0
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	-25	Temperature max	70
Humidity min	5	Humidity max	93



08/12/2022 12:39 - This article is available on our website : www.gsh-med.fr with the following reference SPG445

All characteristics are given for information only and may be changed

2/2