

PVC BLACK TUBE 100 MM

Réf. 520621135



PVC BLACK TUBE 100 MM

Réf. 520621135

FEATURES			
Medical device	No		
Marquage		Barcode	3661040213729
Guarantee (month):	12	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	1,5
Length (cm) :		Length packaging (cm) :	13,0
Height (cm) :		Height packaging (cm) :	1,0
Weight (kg) :		Weight packaging (kg) :	0,018
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:46 - This article is available on our website : www.gsh-med.fr with the following reference 520621135