



TRANSPORT BAG FOR 4VS



Réf. MDSACT01



TRANSPORT BAG FOR 4VS



Réf. **MDSACT01**

FEATURES			
Medical device	No		
Marquage	Barcode	3661040218656	
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	27,0	Width packaging (cm):	
Length (cm) :	42,0	Length packaging (cm) :	
Height (cm) :	15,0	Height packaging (cm) :	
Weight (kg) :	0,500	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:18 - This article is available on our website : www.gsh-med.fr with the following reference MDSACT01

All characteristics are given for information only and may be changed