

Réf. SSARMS01



Réf. SSARMS01

FEATURES			
Medical device	No		
Marquage	Barcode	3661040219875	
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	20,0	Width packaging (cm):	
Length (cm) :	34,0	Length packaging (cm) :	
Height (cm) :	41,4	Height packaging (cm) :	
Weight (kg) :	2,000	Weight packaging (kg) :	2,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:40 - This article is available on our website : www.gsh-med.fr with the following reference SSARMS01

All characteristics are given for information only and may be changed