

Réf. BCASPI142



CE
0123



Réf. BCASPI142

FEATURES			
Medical device	No		
Marquage	CE0123	Barcode	0000000000871
Guarantee (month):	24	Customs Code	90330090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	19,5
Length (cm) :		Length packaging (cm) :	17,0
Height (cm) :		Height packaging (cm) :	7,0
Weight (kg) :		Weight packaging (kg) :	0,272
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:37 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI142

All characteristics are given for information only and may be changed

2/2