

CA-MI NEW ASKIR 30 EKECTRIC SUCTION UNIT POLYCARBONATE JARS 1L

Réf. BCASPI12

CA-MI

GARANTIE
3
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL CLASSE
MEDICAL DEVICE II



CA-MI NEW ASKIR 30 EKECTRIC SUCTION UNIT POLYCARBONATE JARS 1L

Réf. BCASPI12

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	8033622540120
Guarantee (month):	36	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	18,0	Width packaging (cm):	
Length (cm) :	35,0	Length packaging (cm) :	
Height (cm) :	21,0	Height packaging (cm) :	
Weight (kg) :	3,600	Weight packaging (kg) :	3,800
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	1	Humidity max	85



30/03/2022 22:43 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI12