

ELECTRICAL SUCTION UNIT 1000 ML AUTOCLAVABLE POLYCARBONATE JAR

Réf. BCASPI18

CA-MI

GARANTIE
2
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE
II
a



ELECTRICAL SUCTION UNIT 1000 ML AUTOCLAVABLE POLYCARBONATE JAR



Réf. **BCASPI18**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	8054619100370
Guarantee (month):	24	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	35,0	Width packaging (cm):	21,0
Length (cm) :	18,0	Length packaging (cm) :	37,0
Height (cm) :	21,0	Height packaging (cm) :	24,0
Weight (kg) :	3,500	Weight packaging (kg) :	4,600
SPECIFIC STORAGE CONDITION			
Temperature min	-25	Temperature max	70
Humidity min	1	Humidity max	90



30/03/2022 22:43 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI18