

## PORTABLE ELECTRICAL SUCTION UNIT, CX2000ML AUTOCLAVE

Réf. BCASPI30

CA-MI



CE  
0123

DISPOSITIF MEDICAL  
MEDICAL DEVICE

CLASSE  
IIa



# PORTABLE ELECTRICAL SUCTION UNIT, CX2000ML AUTOCLAVE

CA-MI

Réf. BCASPI30

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | Yes       | Classe IIa              |               |
| Marquage                   | CE0123    | Barcode                 | 3661040139166 |
| Guarantee (month):         | 36        | Customs Code            | 90189084      |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               | 30,0      | Width packaging (cm):   |               |
| Length (cm) :              | 32,0      | Length packaging (cm) : |               |
| Height (cm) :              | 99,0      | Height packaging (cm) : |               |
| Weight (kg) :              | 7,890     | Weight packaging (kg) : | 11,600        |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | -20       | Temperature max         | 45            |
| Humidity min               | 0         | Humidity max            | 85            |



30/03/2022 22:43 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCASPI30

All characteristics are given for information only and may be changed

2/2