

CA-MI ASKIR C30 PORTABLE ELECTRICAL SUCTION UNIT, 2X2L WITH WHEELS

Réf. BCASPI34

CA-MI

GARANTIE
2
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE
II



CA-MI ASKIR C30 PORTABLE ELECTRICAL SUCTION UNIT, 2X2L WITH WHEELS



Réf. BCASPI34

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040139203
Guarantee (month):	24	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	32,0	Width packaging (cm):	
Length (cm) :	99,0	Length packaging (cm) :	
Height (cm) :	10,0	Height packaging (cm) :	
Weight (kg) :	6,200	Weight packaging (kg) :	9,800
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 22:44 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI34

All characteristics are given for information only and may be changed