

SUCCION UNIT ASKIR C30, FLOVAC 2X2000ML DISPOSABLE LINERS

Réf. BCASPI37

CA-MI

GARANTIE
3
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE

CLASSE
II
B



SUCCION UNIT ASKIR C30, FLOVAC 2X2000ML DISPOSABLE LINERS

Réf. BCASPI37

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040139234
Guarantee (month):	36	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	54,2	Width packaging (cm):	
Length (cm) :	54,2	Length packaging (cm) :	
Height (cm) :	33,2	Height packaging (cm) :	
Weight (kg) :	9,400	Weight packaging (kg) :	9,400
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 22:44 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI37

All characteristics are given for information only and may be changed