

SUCCION UNIT ASKIRC30,FS2FLOVAC 2X2000ML DISPOSABLE LINERS+FOOTSWITCH

Réf. BCASPI38



CA-MI

GARANTIE
3
ANS
YEARS
GARANTEE

CE
0123

DISPOSITIF MEDICAL CLASSE
MEDICAL DEVICE II

SUCCION UNIT ASKIRC30,FS2FLOVAC 2X2000ML DISPOSABLE LINERS+FOOTSWITCH

CA·MI

Réf. BCASPI38

FEATURES

Medical device Yes Classe IIa

Marquage CE0123

Barcode 3661040139241

Guarantee (month): 36

Customs Code 90189084

Sterile No

INFORMATIONS

Main unit

Packaged unit

Width (cm) : 53,5

Width packaging (cm):

Length (cm) : 55,0

Length packaging (cm) :

Height (cm) : 32,5

Height packaging (cm) :

Weight (kg) :

Weight packaging (kg) : 9,150

SPECIFIC STORAGE CONDITION

Temperature min -20

Temperature max 45

Humidity min 0

Humidity max 85



30/03/2022 22:44 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI38

All characteristics are given for information only and may be changed