

PROFESSIONAL ELECTRICAL SUCTION UNIT, BASIC 2, 2X2L AUTOCLAVABLE JARS


Réf. BCASPI39

CA-MI

GARANTIE
2
ANS
YEARS
GUARANTEED

CE
0123

DISPOSITIF MEDICAL CLASSE II
MEDICAL DEVICE



PROFESSIONAL ELECTRICAL SUCTION UNIT, BASIC 2, 2X2L AUTOCLAVABLE JARS

Réf. **BCASPI39**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040139258
Guarantee (month):	24	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	46,0	Width packaging (cm):	
Length (cm) :	42,0	Length packaging (cm) :	
Height (cm) :	85,0	Height packaging (cm) :	
Weight (kg) :	13,000	Weight packaging (kg) :	19,800
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 22:44 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI39