

## SUCCION UNIT NEW HOSPIVAC 350, 2X2L POLYC JARS, VERSION FULL 2


Réf. BCASPI41

CA-MI

GARANTIE  
**3**  
ANS  
YEARS  
GARANTEE

CE  
0123

DISPOSITIF MEDICAL  
MEDICAL DEVICE **II**



## SUCCION UNIT NEW HOSPIVAC 350, 2X2L POLYC JARS, VERSION FULL 2

CA-MI

Réf. BCASPI41

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | Yes       | Classe IIa              |               |
| Marquage                   | CE0123    | Barcode                 | 3661040139272 |
| Guarantee (month):         | 36        | Customs Code            | 90189084      |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               |           | Width packaging (cm):   | 30,5          |
| Length (cm) :              |           | Length packaging (cm) : | 50,0          |
| Height (cm) :              |           | Height packaging (cm) : | 92,0          |
| Weight (kg) :              |           | Weight packaging (kg) : | 22,000        |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | -20       | Temperature max         | 45            |
| Humidity min               | 1         | Humidity max            | 85            |



30/03/2022 22:44 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCASPI41

All characteristics are given for information only and may be changed