

ROUND SPACER FOR ASKIR 20 AND 30

CA-MI

Réf. BCASPI105



ROUND SPACER FOR ASKIR 20 AND 30

CA-MI

Réf. BCASPI105

FEATURES			
Medical device	No		
Marquage	CE	Barcode	0000000000628
Guarantee (month):	12	Customs Code	90330090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,050
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:27 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI105