

CONICAL CONNECTORS Ø 10-11-12MM FOR ASKIR C30BR

Réf. BCASPI94

CA-MI



CONICAL CONNECTORS Ø 10-11-12MM FOR ASKIR C30BR



Réf. BCASPI94

FEATURES			
Medical device	No		
Marquage		Barcode	3661040169439
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,006
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:27 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI94