

SUCCION UNIT NEW HOSPIVAC 350, FLOVAC 2X2L DISPOSABLE LINERS BASIC 2

Réf. BCASPI48


CA-MI



CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE

CLASSE
II



SUCCION UNIT NEW HOSPIVAC 350, FLOVAC 2X2L DISPOSABLE LINERS BASIC 2

CA-MI

Réf. **BCASPI48**

FEATURES

Medical device Yes Classe IIa

Marquage CE0123

Barcode 3661040139340

Guarantee (month): 24

Customs Code 90189084

Sterile No

INFORMATIONS

Main unit

Packaged unit

Width (cm) : 46,0

Width packaging (cm):

Length (cm) : 42,0

Length packaging (cm) :

Height (cm) : 85,0

Height packaging (cm) :

Weight (kg) : 85,000

Weight packaging (kg) :

SPECIFIC STORAGE CONDITION

Temperature min -20

Temperature max 45

Humidity min 0

Humidity max 85



30/03/2022 22:44 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI48

All characteristics are given for information only and may be changed

2/2

HOLTEX

30 rue Jean de Guiramand
13290 Aix-en-Provence - France

www.gsh-med.fr

Contact :

+33(0)4 42 90 31 31 - info@gsh-med.fr