

S/S BOWL WITHOUT SPOUT, FLAT BOTTOM, DIAM 100 MM, 300 ML, WITH COVER

Réf. AI99101



S/S BOWL WITHOUT SPOUT, FLAT BOTTOM, DIAM 100 MM, 300 ML, WITH COVER



Réf. **AI99101**

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040160573
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	21,5	Width packaging (cm):	21,5
Length (cm) :	38,0	Length packaging (cm) :	38,0
Height (cm) :	4,6	Height packaging (cm) :	4,6
Weight (kg) :	0,096	Weight packaging (kg) :	0,100
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:11 - This article is available on our website : www.gsh-med.fr with the following reference AI99101

All characteristics are given for information only and may be changed