

Réf. BCASPI55

CA-MI

GARANTIE
3
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE IIa



Réf. BCASPI55

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040139418
Guarantee (month):	36	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	50,0
Length (cm) :		Length packaging (cm) :	50,0
Height (cm) :		Height packaging (cm) :	91,0
Weight (kg) :		Weight packaging (kg) :	22,000
SPECIFIC STORAGE CONDITION			
Temperature min	-25	Temperature max	70
Humidity min	1	Humidity max	93



30/03/2022 22:45 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI55

All characteristics are given for information only and may be changed

2/2