

ELECTRICAL SUCTION UNIT 2X5000ML, NEW HOSPIVAC 400, VERSION FULL

Réf. BCASPI59

CA-MI



CE
0123



ELECTRICAL SUCTION UNIT 2X5000ML, NEW HOSPIVAC 400, VERSION FULL

CA-MI

Réf. **BCASPI59**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040139456
Guarantee (month):	36	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 22:45 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI59

All characteristics are given for information only and may be changed

2 / 2