

ANTIBACTERIAL & HYDROPHOBIC FILTER FOR ASKIR C30BR

Réf. BCASPI97



CA-MI



ANTIBACTERIAL & HYDROPHOBIC FILTER FOR ASKIR C30BR

CA-MI

Réf. **BCASPI97**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	8054619103869
Guarantee (month):	12	Customs Code	90192000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	16,0	Width packaging (cm):	
Length (cm) :	25,0	Length packaging (cm) :	
Height (cm) :	5,5	Height packaging (cm) :	
Weight (kg) :	0,050	Weight packaging (kg) :	0,050
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:01 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI97

All characteristics are given for information only and may be changed

2/2