

DISPOSABLE FLOVAC® LINER, 1000 ML FOR ASPIRET, ASKIR 20, 30

Réf. BCASPI99

CA-MI

GARANTIE
1
AN
TOUT
RISQUE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE
I_m



DISPOSABLE FLOVAC® LINER, 1000 ML FOR ASPIRET, ASKIR 20, 30

CA-MI

Réf. BCASPI99

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0123	Barcode	3661040151137
Guarantee (month):	12	Customs Code	90189084
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		0,200	Weight packaging (kg) : 0,200
SPECIFIC STORAGE CONDITION			
Temperature min		-20	Temperature max 45
Humidity min		0	Humidity max 85



30/03/2022 23:01 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI99

All characteristics are given for information only and may be changed

2/2