

SIMS SCISSORS B/S STRAIGHT 20 CM

Réf. IC23220



SIMS SCISSORS B/S STRAIGHT 20 CM



Réf. IC23220

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040010809
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	5,3	Width packaging (cm):	
Length (cm) :	20,3	Length packaging (cm) :	
Height (cm) :	0,7	Height packaging (cm) :	
Weight (kg) :	0,078	Weight packaging (kg) :	0,078
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:06 - This article is available on our website : www.gsh-med.fr with the following reference IC23220

All characteristics are given for information only and may be changed