

CA-MI NEW ASKIR 30 PROXIMITY ELECTR SUCTION UNIT, 1L JAR IN MAKROLON

Réf. BCASPI400

CA-MI

GARANTIE
2
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE
II



CA-MI NEW ASKIR 30 PROXIMITY ELECTR SUCTION UNIT, 1L JAR IN MAKROLON



Réf. **BCASPI400**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	8054619100400
Guarantee (month):	24	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	35,0	Width packaging (cm):	
Length (cm) :	21,0	Length packaging (cm) :	
Height (cm) :	18,0	Height packaging (cm) :	
Weight (kg) :	3,600	Weight packaging (kg) :	4,250
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:33 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI400