


"O" RING FOR JAR 1000 AND 2000 ML

Réf. BCASPI87

CA-MI



"O" RING FOR JAR 1000 AND 2000 ML

CA-MI

Réf. BCASPI87

FEATURES			
Medical device	No		
Marquage		Barcode	3661040169361
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,200	Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



15/11/2022 12:44 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI87

All characteristics are given for information only and may be changed