

"O" RING FOR JAR 4000 ML

Réf. BCASPI89



"O" RING FOR JAR 4000 ML



Réf. BCASPI89

FEATURES			
Medical device	No		
Marquage		Barcode	3661040169385
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,200	Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:34 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI89