

FOOT SWITCH FOR ASKIR C30/BR AND HOSPIVAC 350/400

CA-MI

Réf. BCASPI90



FOOT SWITCH FOR ASKIR C30/BR AND HOSPIVAC 350/400

CA-MI

Réf. **BCASPI90**

FEATURES			
Medical device	No		
Marquage		Barcode	3661040169392
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		0,200	
		Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min		-20	
		Temperature max	45
Humidity min		0	
		Humidity max	85



30/03/2022 23:34 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI90

All characteristics are given for information only and may be changed

2/2