

## WEITLANER RETRACTOR 10.5 CM

Réf. IE31110



## WEITLANER RETRACTOR 10.5 CM



Réf. IE31110

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040017235
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,3	Width packaging (cm):	
Length (cm) :	11,2	Length packaging (cm) :	
Height (cm) :	2,0	Height packaging (cm) :	
Weight (kg) :	0,052	Weight packaging (kg) :	0,052
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



15/11/2022 21:11 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IE31110

All characteristics are given for information only and may be changed