

WEITLANER RETRACTOR 24 CM (5X6 PRONGS BLUNT)

Réf. IE31124



WEITLANER RETRACTOR 24 CM (5X6 PRONGS BLUNT)



Réf. IE31124

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe IIa | |
| Marquage | CE | Barcode | 3661040017273 |
| Guarantee (month): | 24 | Customs Code | 9018.9090 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 9,0 | Width packaging (cm): | |
| Length (cm) : | 24,4 | Length packaging (cm) : | |
| Height (cm) : | 3,2 | Height packaging (cm) : | |
| Weight (kg) : | 0,188 | Weight packaging (kg) : | 0,188 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



30/03/2022 16:38 - This article is available on our website : www.gsh-med.fr with the following reference IE31124

All characteristics are given for information only and may be changed