

## STAINLESS STEEL INSTRUMENT TRAY 40X30X4 CM

Réf. AI40304



# STAINLESS STEEL INSTRUMENT TRAY 40X30X4 CM



Réf. AI40304

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040067919
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	30,0	Width packaging (cm):	35,5
Length (cm) :	40,0	Length packaging (cm) :	46,0
Height (cm) :	4,0	Height packaging (cm) :	6,5
Weight (kg) :	0,900	Weight packaging (kg) :	1,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:49 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference AI40304

All characteristics are given for information only and may be changed