

## S/S BOX 50 X 20 X 12 CM

Réf. AI50201



## S/S BOX 50 X 20 X 12 CM



Réf. AI50201

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040019772
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	20,0	Width packaging (cm):	50,0
Length (cm) :	50,0	Length packaging (cm) :	20,0
Height (cm) :	12,0	Height packaging (cm) :	12,0
Weight (kg) :	2,263	Weight packaging (kg) :	2,263
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:49 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference AI50201

All characteristics are given for information only and may be changed