

LECOUTEUR TONGUE DEPRESSOR

Réf. IA02420



LECOUTEUR TONGUE DEPRESSOR



Réf. IA02420

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040006000
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	2,7	Width packaging (cm):	
Length (cm) :	16,0	Length packaging (cm) :	
Height (cm) :	7,5	Height packaging (cm) :	
Weight (kg) :	0,051	Weight packaging (kg) :	0,051
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:41 - This article is available on our website : www.gsh-med.fr with the following reference IA02420

All characteristics are given for information only and may be changed