


REDON INTRODUCING NEEDLE STRAIGHT (DIAM. 3.5 MM), 15CM

Réf. IA08410



REDON INTRODUCING NEEDLE STRAIGHT (DIAM. 3.5 MM), 15CM



Réf. IA08410

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040006413
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	11,0	Width packaging (cm):	
Length (cm) :	15,8	Length packaging (cm) :	
Height (cm) :	0,3	Height packaging (cm) :	
Weight (kg) :	0,018	Weight packaging (kg) :	0,018
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:43 - This article is available on our website : www.gsh-med.fr with the following reference IA08410

All characteristics are given for information only and may be changed