

HYPOPHYSECTOMY CANNULA STRAIGHT (SET OF 4 Ø1MM)

Réf. IC41110



HYPOPHYSECTOMY CANNULA STRAIGHT (SET OF 4 Ø1MM)



Réf. IC41110

FEATURES			
Medical device	No		
Marquage		Barcode	3661040110592
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,032	Weight packaging (kg) :	0,032
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:12 - This article is available on our website : www.gsh-med.fr with the following reference IC41110

All characteristics are given for information only and may be changed