

LUER ST STEEL TRACHEAL TUBE 2 PCS N°3 DIAM. 7 MM

Réf. IC43203



LUER ST STEEL TRACHEAL TUBE 2 PCS N°3 DIAM. 7 MM



Réf. IC43203

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040011653
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	4,4	Width packaging (cm):	
Length (cm) :	7,3	Length packaging (cm) :	
Height (cm) :	2,4	Height packaging (cm) :	
Weight (kg) :	0,018	Weight packaging (kg) :	0,018
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



25/05/2022 12:33 - This article is available on our website : www.gsh-med.fr with the following reference IC43203

All characteristics are given for information only and may be changed